

# Advanced Characterizations of Gas-Polymer Interactions: Pioneering Carbon Capture and Storage Research into the ECCSELLENT Project


Virginia Signorini<sup>1</sup>, Gianfranco Burzotta<sup>1</sup>, Roberta Di Carlo<sup>1</sup>, Zahra Maghazeh<sup>1</sup>, Marco Giacinti Baschetti<sup>1</sup>, Matteo Minelli<sup>1</sup>,

<sup>1</sup>Departement of Civil, Chemical, Environmental and Materials engineering, University of Bologna, via Terracini 28, Bologna

Carbon Capture and Storage (CCS) is a crucial technology aimed at reducing greenhouse gas emissions by capturing CO<sub>2</sub> directly from industrial sources and power plants [1,2]. Within the **ECCSELLENT** project (ECCSEL framework), the University of Bologna is committed to improve the Italian research infrastructure in the **CCUS** field (CO<sub>2</sub> Capture, Utilization, Transport, and Storage), with particular focus on CO<sub>2</sub> capture by membrane and on CO<sub>2</sub> transportation. Equipped with state-of-the-art instrumentation, **MEMLAB** and **Transport LAB** enables comprehensive characterization of polymers, evaluating their transport properties—including permeability, solubility and diffusivity—under a broad range of operating conditions.

## MEMLAB

Polymeric membrane morphology and nanostructure characterization for CO<sub>2</sub> capture capacity, separation and purification performances, in both pre- and post-combustion scenarios.




### FT-IR




- Interferometer, laser and IR source for streamline data analysis.
- Diamond crystal ATR accessory.
- Temperature controlled plate (273K to 573K).
- Pressure dome from 1 to 70bar

### QCM-D



- Real-time frequency and energy dissipation
- Mass adsorption
- Viscoelastic properties
- Gas and Liquid environments


### ATOMIC FORCE MICROSCOPY



- High-resolution surface characterization
- Static and dynamic imaging modes
- Air and liquid environments
- Force-distance and tip current-voltage mapping.

## TRANSPORT LAB

Characterization of gas and liquid transport through polymeric materials to determine polymer-penetrant interaction, even in cryogenic and supercritical conditions, for transport applications.



### HIGH-PRESSURE SORPTION ANALYZER



- Automated gas sorption analyser for volumetric measurements.
- Recording adsorption and desorption isotherms
- Temperature: 20K to 773K
- Pressure: 1 bar to 200 bar

### DIELECTRIC SPECTROMETER




- Determination of polymeric chain mobility and relaxation.
- CO<sub>2</sub> effects on material's conductivity and permittivity.
- Temperature: 113K to 673K

### HIGH PRESSURE PERMEOMETER





- Fixed volume and variable pressure permeometer
- Determination of Permeability (P) and Diffusion Coefficient (D)
- Temperature: 260K to 338K
- Pressure up to 150 bar.




## THE RESULTS

CO<sub>2</sub> desorption in Natural Rubber determined through FT-IR Spectrometer




Investigation of Matrimid dielectric properties upon CO<sub>2</sub> exposure at 40 bar



Surface, Phase-segregation, roughness and morphological characterization of PS-PEO Block-copolymers

CO<sub>2</sub> uptake isotherms in PIM-1 as a function of gas activity, at various cryogenic temperatures



## THE OUTCOMES


Thermodynamic EoS (Lattice Fluid / Non-Equilibrium Lattice Fluid (NELF)) are used to describe the effect of CO<sub>2</sub> on polymers across a range of temperatures and pressures relevant to CO<sub>2</sub> transport and capture applications.

$$\mu_i^{pol}(T, p, \omega_i, p_{pol}) = \mu_i^{gas}(T, p, \omega_i)$$


The **Standard Transport Model (STM)** integrates solubility and diffusion coefficient to describe the gas permeability behavior as a function of pressure [5].

$$P_i = \frac{1}{(p_i^u - p_i^d)} \int_{p_i^d}^{p_i^u} L_0 e^{\beta \Omega_i} S_i z_i dp_i$$

$L_0$ : mobility coefficient  
 $\beta$ : plasticization factor  
 $S_i$ : solubility coefficient  
 $z_i$ : penetrant compressibility  
 $dp_i$ : pressure difference



This modeling framework enables a deeper understanding of polymer-penetrant interactions and the effects of dense-phase CO<sub>2</sub> on polymer materials, while predicting polymer behavior across a range of temperatures and pressures. As a result, it supports the rational design and optimization of materials for CO<sub>2</sub> capture and transport applications.



## REFERENCE

- [1] D.Y.C. Leung, G. Caramanna, M.M. Maroto-Valer, Renew. Sustain. Energy Rev. (2014)
- [2] L. Ansaloni, B. Alcock, T.A. Peters, nt. J. Greenh. Gas Control. 94 (2020) 102930.
- [3] F. Doghieri, G.C. Sarti, Nonequilibrium lattice fluids, Macromolecules. 29 (1996) 7885–7896.
- [4] I.C. Sanchez, R.H. Lacombe, Macromolecules 11 (1978) 1145–1156.
- [5] E. Ricci, M.G. De Angelis, M. Minelli, Chem. Eng. J. 435 (2022) 135013.



## ACKNOWLEDGMENTS

National Recovery and Resilience Plan, Mission 4 Component 2 Investment 3.1 "Fund for the realisation of an integrated system of research and innovation infrastructures". Call for tender 3264 (28 Dec. 2021) of the Italian Ministry of University and Research funded by the European Union—NextGenerationEU—PNRR IR0000020, Concession Decree 244 (8 Aug. 2022) adopted by the Italian Ministry of University and Research, CUP F53C22000560006, ECCSELLENT—Development of ECCSEL-R.I. Italian facilities: user access, services and long-term sustainability.