

Stefan Helfert,^{1,6} Chrissie I.M. Baltzaki^{1,2,6}, Tommaso Zandrini,^{3,6} Andreas Rohatschek,^{2,4,6} Manuel Rufin,^{4,6} Peter Machata,⁵ Anna Zahoranová,^{1,2,6*} Orestis G. Andriotis,^{4,6} Philipp J. Thurner,^{4,6} Aleksandr Ovsianikov,^{3,6} Robert Liska,^{1,6} Stefan Baudis^{1,2,6*}

¹ Institute of Applied Synthetic Chemistry, Technische Universität Wien, Vienna, Austria, ² Christian Doppler Laboratory for Advanced Polymers for Biomaterials and 3D Printing, Technische Universität Wien, Vienna, Austria, ³ Institute of Materials Science and Technology, Technische Universität Wien, Vienna, Austria

⁴ Institute of Lightweight Design and Structural Biomechanics, Technische Universität Wien, Vienna, Austria, ⁵ Department of Composite Materials, Polymer Institute of the Slovak Academy of Sciences, Bratislava, Slovakia, ⁶ Austrian Cluster for Tissue Regeneration, Vienna, Austria



Contact: chrysi.baltzaki@tuwien.ac.at, www.baudislab.com




Austrian Cluster for Tissue Regeneration

Introduction


Polymer brushes can provide tunable surface functionality –vital for sensors, implants and tissue scaffolds. Challenge: Patterning polymer brushes in complex geometries limited by current methods.

Solution: 2-Photon absorption to initiate RAFT polymerization directly on surfaces



Methods


- RAFT agent: **CDTPA**
- Monomer: Biobiocompatible & hydrophilic **NAM**
- Substrate: Si wafers modified with CDTPA
- Conditions: 420:1 NAM:CDTPA, blue light



2Photon RAFT (2PRAFT)

Transition from RAFT to 2PRAFT:



- Change of photoinitiator to M2CMK
- Formulation: NAM:CDTPA:FloAc:M2CMK=420:1:1:1
- Substrate: Microscopy Coverslips
- Optimal Conditions:**
 - Fluorescence @ 100 mm/s
 - Best performance @ 80 mW



Multichemistry & Confined Printing

Multicolour: FloAc & Nile Blue acrylamide (**NBAAm**)

Formulations:
NAM:CDTPA:FloAc:M2CMK=420:1:1:1
NAM:CDTPA:NBAAm:M2CMK=420:1:1:1



- ✓ Dual-side patterning
- ✓ No cross-talk between layers
- ✓ Fine resolution up to 5 µm



Conclusions

- ✓ 2PRAFT method was developed for the first time, and used for polymer brush patterning.
- ✓ Well-defined brushes with 10 nm height and structures with edges down to 5 µm were produced.
- ✓ Created multicolour images & print different patterns on adjacent surfaces.
- ✓ Applications: Controlled cell adhesion, microfluidic devices.

REFERENCES & ACKNOWLEDGEMENTS

[1] Helfert, S.; Zandrini, T.; Rohatschek, A.; Rufin, M.; Machata, P.; Zahoranová, A.; Andriotis, O.G.; Thurner, P.J.; Ovsianikov, A.; Liska, R.; Baudis, S.; Micropatterning of Confined Surfaces with Polymer Brushes by Two-Photon-Initiated Reversible Addition–Fragmentation Chain-Transfer Polymerization. *Small Science*, 5(1), p.2400263, DOI: 10.1002/smssc.202400263.

Funding by the Christian Doppler Research Association within the framework of a Christian Doppler Laboratory for “Advanced Polymers for Biomaterials and 3D Printing” and the financial support by the Austrian Federal Ministry of Labor and Economy Affairs and the National foundation for Research, Technology and Development are gratefully acknowledged.

More info:

